

機器レイアウト



ITEM	DESCRIPTION	DIMENSIONS LxWxH (mm)	WEIGHT (kg)
1	LCフォーム	H2234	750
2	LPフォーム	-	797
3	患者撮影ロングテーブル	727x3336x780	590
4	天井走行ルール	4724x64x84	5.1kg/m
5	LC用カーゴバルブレイアール	-	-
6	TVモニター用支持架	2900x635x162	66
7	LargeDisplay ※	144x1580x900	356
8	X線防護板用ルール ※	300x2500x70	69
9	C1キヤビネット (F)	654x900x1854	404
10	C1キヤビネット (L)	654x900x1854	320
11	C2キヤビネット	654x900x1854	283
12	X線管チラーユニット(Coolix4100)	610x555x1200	120
13	Autotransformer (Coolix4100)	304x370x340	30
14	ダイナタタチラーユニット	350x200x340	15
15	UPS 3KVA	412x214x346	35
16	シスナデルPDB (GEタイプ)	300x800x1300	225
17	J-PDB	400x900x1550	245
18	ドキュメント用ロッカー	-	-
19	分電盤 (お客様準備品)	600x300x800	42
20	専用デスク	700x1800x700	-
21	AdvantageWorkStation ※	-	37
22	オペレーターコンソール	-	-
23	リフレクソスモニター	-	-
24	LargeDisplayCabinet (LD付属) ※	780x550x1250	115
25	油圧台車	900x450x950	-
26	CAPフロントム	400x600x470	-
27	3Dフロントム	580x800x460	-
28	3KVA UPS (LargeDisplay Monitor)	214x421x346	34.5

天井仕上げ面高さ(撮影室) Fix: 2845mm ± 5mm

社会医療法人壮幸会 行田総合病院 御中
JAPAN

GE Healthcare
Satoshi Mori
Satoshi.mori@ge.com

IGS 630
PRELIMINARY STUDY

ピット、棚、医ガス、フィルター他追記
図面修正
新規作成
A 12/Jan/2018
DATE

設置場所要検討
26 27 18

GE does not take responsibility for any damages resulting from changes or drawings made by others. Errors may occur by not referring to the complete set of final issue drawing. GE cannot accept responsibility for any damages due to the partial use of GE final issue drawings, however caused. Do not scale from printed pdf files. GE accepts no responsibility or liability for defective work due to scaling from these drawings.

A mandatory component of this drawing set is the GE Healthcare Pre installation manual. Failure to reference the Pre installation manual will result in incomplete documentation required for site design and preparation. Pre installation documents for GE Healthcare products can be accessed on the web at: www.gehealthcare.com/sipanning

Drawn by	Verified by	S.O. (GON)	Concession	PIM Manual	Rev
Kaneko	Tanaka	-	-	5694389-1-1EN	2

Format	Scale	File Name	Date	Rev	Sheet
A3	1:60	IGS-A043996-CDWG	02/Feb/2018	C	01/02